

## Goat Anti-Affibody<sup>®</sup> IgG

### BACKGROUND

The goat Anti-Affibody<sup>®</sup> IgG reagent recognizes and binds to Affibody<sup>®</sup> molecules. The antibody has been purified by affinity chromatography and is available as an unconjugated reagent. It can be used as a second step reagent together with a labeled anti-goat antibody in applications where an Affibody<sup>®</sup> molecule is used as a primary reagent, such as immunofluorescence, immunohistochemistry and flow cytometry applications.

Affibody<sup>®</sup> products are supplied for research use only and may not be used in a manufacturing process or in manufactured products or in human subjects. The Affibody<sup>®</sup> products shall in no event be used in medical or clinical applications, unless otherwise expressly stated by Affibody. The purchaser of Affibody<sup>®</sup> products is solely liable to ensure compliance with any regulatory requirements related to the purchaser's use of Affibody<sup>®</sup> products.

### PRODUCT INFORMATION

**Product name:** Goat Anti-Affibody<sup>®</sup> IgG

**Catalog number:**

500 µl: 20.1000.01.0005

**Source:** Polyclonal Anti-Affibody<sup>®</sup> antibodies were raised in goats by repetitive immunization using Affibody<sup>®</sup> molecules as immunogens. Specific anti-Affibody<sup>®</sup> IgG was purified from serum by antigen specific affinity chromatography.

**Specificity:** Goat Anti-Affibody<sup>®</sup> IgG binds to Affibody<sup>®</sup> molecules.

**Applications:** This secondary antibody can be used to detect Affibody<sup>®</sup> molecules.

**Form:** In buffer solution, PBS pH 7.4, in 50% glycerol.

**Storage:** -20°C is recommended.

**Shipping:** At ambient temperature.

**Stability:** There is no decrease in performance of the product after 4 weeks at room temperature.

**Product support:** [www.affibody.com/shop](http://www.affibody.com/shop)  
Affibody AB, P.O. Box 20137, SE-161 02 Bromma, Sweden. Phone: +46-8-598 838 00, Fax: +46-8-598 838 01. E-mail: [biotechnology@affibody.com](mailto:biotechnology@affibody.com).

### LIMITATIONS

Warranty: Affibody<sup>®</sup> products are warranted to meet stated product specifications and to conform to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to one year from date of sales for products used, handled and stored according to Affibody AB's instructions. Affibody AB's sole liability is limited to replacement of the product or refund of the purchase price. Affibody<sup>®</sup> products are supplied for research use only. They are not intended for medicinal, diagnostic or therapeutic use. Affibody<sup>®</sup> products may not be resold, modified for resale or used to manufacture commercial products without prior written approval from Affibody AB.

AFFIBODY AB, PO Box 20137, SE-161 02 Bromma, Sweden. Phone: +46-8-598 838 00, Fax: +46-8-598 838 01, E-mail: [biotechnology@affibody.com](mailto:biotechnology@affibody.com), Web: [www.affibody.com/shop](http://www.affibody.com/shop).  
Rev 070416