

Anti-EGFR Affibody[®] Molecule, Agarose Immobilized

BACKGROUND

The Epidermal Growth Factor Receptor (EGFR, ErbB) is a cell surface glycoprotein of approximately 135 kD (unglycosylated). There are several alternatively spliced forms of EGFR including a secreted soluble form. The EGFR belongs to the family tyrosine kinase receptors which are characterized by an extracellular ligand-binding domain, a single transmembrane and an intracellular domain responsible for transducing the signal. The receptor dimerizes upon binding to EGFR and the transduced signal stimulates cell growth and differentiation. The EGFR is over expressed or mutated in many types of cancers and the receptor as well as the down stream signaling molecules are subjects to several cancer therapeutic interventions.

The Anti-EGFR Affibody[®] molecule is a specific affinity ligand selected against the extracellular domain of EGFR. The agarose immobilized Anti-EGFR Affibody[®] molecule can be used for immunoprecipitation of EGFR from complex mixes of proteins.

PRODUCT INFORMATION

Product name: Anti-EGFR Affibody[®] molecule, agarose immobilized

Catalog number:
200 µl: 10.1246.04.0002
1 ml: 10.1246.04.0010

Source: Recombinant protein produced in *E. coli*.

Specificity: Anti-EGFR Affibody[®] molecule binds to human EGFR. Cross reactivity with other species has not been tested.

MW: 13.7 kDa

Purity: >98% as determined by SDS-PAGE and RP-HPLC analysis.

Tested applications: Immunoprecipitation, Affinity Chromatography.

Form: As 50% slurry in buffer solution (20 mM Sodium Phosphate, 150 mM Sodium Chloride, 0.02%(w/v) Sodium Azide, pH 7.5).

Storage: At +4°C. Avoid freezing.

Shipping: At ambient temperature.

Stability: There is no decrease in performance of the agarose immobilized Affibody[®] molecule after storage for 2 weeks at +37°C.

Product support: www.affibody.com/shop
Affibody AB, P.O. Box 20137, SE-161 02 Bromma, Sweden. Phone: +46-8-598 838 00, Fax: +46-8-598 838 01. E-mail: biotechnology@affibody.com.

LIMITATIONS

Warranty: Affibody[®] products are warranted to meet stated product specifications and to confirm to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to one year from date of sales for products used, handled and stored according to Affibody AB's instructions. Affibody AB's sole liability is limited to replacement of the product or refund of the purchase price. Affibody[®] products are supplied for research use only. They are not intended for medicinal, diagnostic or therapeutic use. Affibody[®] products may not be resold, modified for resale or used to manufacture commercial products without prior written approval from Affibody AB.

AFFIBODY AB, PO Box 20137, SE-161 02 Bromma, Sweden. Phone: +46-8-598 838 00, Fax: +46-8-598 838 01, E-mail: biotechnology@affibody.com, Web: www.affibody.com/shop.
Rev 061027